
Exercise Instructions for Galaxy Tutorial


- 1. Get data from EBI SRA exercise
- a. Click on Get data on the left side panel
- b. Click on EBI SRA
- c. Enter SRR1016067 on the search box and
- d. Click search
- e. Click on Fastq files (galaxy) File1
- f. repeat steps a-d and select c click on Fastq files (galaxy) File2




g. check if you were redirected back to maxwell.abdn.ac.uk/galaxy and not usegalaxy.org, or try again




- h. click on the pencil icon for result 1 in history
- i. change the name to SRR1016067_1.fastq (or any name without gz extension)
- j. select "SacCer2" in the Database/Build" dropdown and save
- k. click Save




- l. click on the pencil icon for result 2 in history
- m. change the name to SRR1016067_2.fastq (or any name without gz extension)
- n. select "SacCer2" in the Database/Build" dropdown and save
- o. click Save




- 2. QA using FastQC exercise
- a. click on NGS:QC and manipulation
- b. click on FASTQC: Reads QC
- c: Select result 1
- d. Click execute



e. click on the eye icon to visualize the results f. scroll down to per base sequence




- 3. Filtering using trimgalore! exercise
- a. click on Trim Galore in the tools panel
- b click on trimgalore sub menu
- c. select paired-end
- d. select result 1
- e. select result 2
- f. check that all parameters are adequate
- g. click on execute



h. click on NGS: Qc and manipulation

i. click on FASTQC: Reads QC

j: Select result 4 k. Click execute




l. click on the eye icon to visualize the results m. scroll down to per base sequence




4a. Upload data using ftp exercise

- a. click on Get Data
- b. click on Upload file
- c. select File Format fasta
- d. paste into URL/File ftp://ftp.ensembl.org/pub/release-
- $59/fasta/saccharomyces_cerevisiae/dna/Saccharomyces_cerevisiae. EF 2.59. dna. toplevel. fa. gz$
- e. select Genome SacCer2
- f. click Execute



4b. Modify datatype


- a. click on the pencil icon for Trim Galore trimmed reads pair 1
- b. select New Type: fastqsanger
- c. Click Save
- d. click on the pencil icon for Trim Galore trimmed reads pair 2
- e. select New Type: fastqsanger
- f. Click Save



4c. Bowtie alignment


- a. Click on Bowtie2 category
- b. Click on Bowtie2
- c. Select paired-end
- d. Select Trim Galore trimmed reads pair 1
- e. Select Trim Galore trimmed reads pair 2
- f. Select "use one form history"
- g. Select result 7 (ftp://ftp.ensembl.org/pub/release-
- 59/fasta/saccharomyces_cerevisiae/dna/Saccharomyces_cerevisiae.EF2.59.dna.toplevel.fa)

h. Click Execute




4d. Visualizing results in UCSC


- a. Click on result 8 link (on the underlined result name "Bowtie 2 on data...")
- b. Click on display at UCSC Main




- d. Click "refresh"
- e. Click Zoom 1.5x




f. use the buttons to navigate g. or right click on the browser to navigate




- 5a. Sharing history and creating workflow exercises (requires Galaxy Account)
- a. click on User -> Login
- b. Enter your meail address
- c. enter your password
- d. Click "Login"



- e. Click on saved histories
- f. Click on rename



- g. Enter a name for the History
- h. Click rename history




i. click on share or publish

j. click Make accessible and publish



- 5b. Create workflow exercise (requires Galaxy Account)
- a. Select Extract Workflow
- b. Uncheck "treat as input dataset" 3x
- c. Click on "Create Workflow"




- 6. Converting BAM to UCSC bedGraph using a public workflow
- a. click on Shared Data -> Published Workflows



b. Click on Ensembl BAM to UCSC bedGraph






d. Click on imported: ensembl BAM to UCSC bedGraph




- e. Select result 8 (Bowtie 2 on data...")
- f. Click on run Workflow



- g. click on the pencil icon
- h. Select Datatype tab
- i. Select New Type: bedGraph and press Save



- i. Click on the result link
- k. Click on display at UCSC Main



- l. click on fullm. click refresh
- Genomes Genome Browser **Tools** Mirrors Downloads My Data View Help About Us UCSC Genome Browser on S. cerevisiae June 2008 (SGD/sacCer2) Assembly move <<< <> >> >>> zoom in 1.5x 3x 10x base zoom out 1.5x 3x 10x 100x chrl:12,510-35,633 23,124 bp. enter position or search terms l saccer2 15,000 User Supplied Track User Track воні 🚺 Y6L 864C-8 Other Features from Saccha YALWdeltai 🛐 FLO9 🐼 Human Proteins Mapped by Chained tBLASTn MUC17 ⟨\ MUC17 ⟨\ MK095698 ⟨\ NUFIP2 Non-S. cerevisiae RefSeq Genes Other RefSea Publications: Sequences in Scientific Articles
 S. cerevisiae mRNAs from GenBank S. cerevisiae ESTs That Have Been Spliced Spliced ESTs Eran Segal Regulatory Module Reg. Module Conservation sacPar sacMik | | | | sacKud | | | | sacBay | | | | sacCas sacKlu Click on a feature for details. Click or drag in the base position track move start move end to zoom in. Click side bars for track options. Drag side bars or labels < 2.0 up or down to reorder tracks. Drag tracks left or right to new position. default order hide all manage custom tracks track hubs configure reverse resize refresh [track search] [default tracks] Use drop-down controls below and press refresh to alter tracks displayed. collapse all expand all Tracks with lots of items will automatically be displayed in more compact modes **Custom Tracks** refresh hide dense

n. click on the bowser area to scroll along the genome

